



For the effective treatment of mild to moderate acne.

Omnilux clear-U™ hand-held has been developed from the Omnilux Blue medical device that has been used by dermatologists worldwide for the successful treatment of acne. Omnilux clear-U™ incorporates blue and red LED light with the flick of the switch. Omnilux clear-U™ uses light technology in the form of light emitting diodes to:

- Destroy the bacteria responsible for acne
- Stimulate healthy tissue growth
- Decrease inflammation and pore size
- Limit oil production in the future

Clear-U benefits anyone suffering with mild to moderate inflammatory acne by killing the key bacteria responsible for acne - Propionibacterium acnes (P.acnes). This bacteria produces natural chemicals called porphyrins, which are sensitive to light at specific wavelengths. When your skin is flooded with blue and red light, the bacteria are killed and inflammation is decreased. Acne fades and skin heals. It's as simple as that!

When?

The Omnilux clear-U™ should be used at the first sign of an acne breakout. The treatment program consists of 8 treatment sessions spread over 4 weeks, 2 treatments per week. Use the blue light on the first treatment session; wait 2 days, then use the red light. Repeat this for the next 3 weeks.

What results can be expected from a series of treatments?

Best results will be seen between 4-8 weeks after treatment. The light stimulates natural processes that occur after the treatment has stopped. Studies have shown complete clearance of lesions, however, on average you should expect 7 out of 10 acne lesions to clear.

How many treatments are needed?

A recommended course of treatments is 2x20 minute treatments per week for 4 weeks of alternating blue and red light. However, because the light will not damage your skin, there is no limit on the number of treatments you can have.

Effects last for six months. In some cases, follow-up treatments (typically 2 to 4 treatments) may be necessary.

clear-U is painless and has only produced mild redness that recedes within one to two days.



UNIT SIZE
H = 12 cm W = 6 cm
D = 4 cm at maximum

Omnilux clear-U™ is not a substitute for good skin hygiene or a doctor's care. We suggest that the light therapy is used in conjunction with a preparation containing salicylic acid.

Omnilux clear-U™ Recommended Treatment Protocol
Follow these simple instructions for clearer skin:

- Cleanse the face using an oil-free, fragrance free cleansing milk or wash
- Use Omnilux clear-U™ as directed
- Apply an appropriate moisturiser to the face
- Follow with an SPF 30+ sunscreen during daylight hours
- Wait at least 48 hours before your next clear-U light treatment

PRODUCT:

Omnilux clear-U™ hand-held device RRP \$690

Australian Stockists: 1800 808 993